


The Torch™ System UV Data Sheet

ClorDiSys Solutions' Torch™ UV Disinfection System is designed to be the most efficient, effective, and affordable UV-C decontamination system for rooms available. The Torch is a portable, high efficiency device that is designed for reliable daily use to reduce organisms at your facility. On the outside, stainless steel construction and large diameter casters provide a long lasting device that is easy to maneuver. On the inside, the highest quality ballasts and bulbs, coupled with an industrial PLC controller ensure years of trouble free use. To keep the cost down, versions are offered with and without special features so no one needs to pay for more than what they want. The Torch system itself offers the lowest price per UV-C watt output available. Treatment costs under a dollar and the low cost of lamps further enhances the affordability.

UV output was designed to obtain greater than 99% reduction of typical viruses and bacteria in a 1-minute timeframe and on spores like C. diff in a 5-minute timeframe within an 8 ft distance. The Torch produces a UVC intensity of approximately 12 mJ/cm² per minute (200 μw/cm²) at an 8 feet distance.



Distance	UV-C Intensity	Dosage per minute
4 ft	378 μw/cm ²	22.68 mJ/cm ²
8 ft	200 μw/cm ²	12 mJ/cm ²
10 ft	128 μw/cm ²	7.68 mJ/cm ²

Background:

- On average, 5% of hospital patients develop an HAI, and 10% of ICU patients develop HAIs (Grohskopf 2002).
- 1.7 Million Americans contract an HAI every year. 99,000 of these patients die from the complications of an HAI (Srinivasan 2009).
- On average, hospital stays with infections due to medical care were 19.2 days longer and the cost was nearly \$43,000 greater than stays without infections (Lucado 2010).
- CAUTIs, the most common HAI, account for 33% of all HAIs. CLABSIs and VAPs each account for 14%, meaning >60% are device-associated (Cass 2013)
- 33% of Operating Rooms responded as to having an infection or outbreak in the last six months (ICT 2013).



- 41% of patient rooms had at least one surface contaminated with MRSA and/or C. difficile (Faires 2013).
- Up to 60% of hospital uniforms are colonized with potentially pathogenic bacteria (Wiener-Well 2011).
- For healthcare workers entering a room containing a patient with MRSA infection, the bacteria would be found on the healthcare worker's clothes approximately 70 percent of the time, even if the healthcare worker did not touch the patient (Pyrek 2012).

Features:

Effective:

- Multiple Torch systems can be utilized to treat larger areas or areas with complex shapes to get optimal coverage. A typical use consists of one unit in the main part of the room and (for hospital patient rooms) another in the bathroom. It must be understood that all UV light works by "line-of-site" meaning that the surfaces need to be illuminated by the light in order to achieve effective kill. Shadowed areas are not affected by the UV Light. Basically two Torches might have the same UV-C output as competing systems but get better disinfection because overall coverage is better with multiple units. The cost of 2 Torch units is less than the cost of the larger competing systems. The Torch was designed to optimize its UV-C output vs. power usage and cost, so that multiple Torches can be placed into a room to get more effective kill than one competitor unit with double the output and more than twice the cost.
- Eight high-output UV-C bulbs are utilized to get optimal intensity -- balanced with power usage -- for efficient kill.
- The center of the Torch is open so that each of the 8 UV-C bulbs can radiate its light 360 degrees.
- The UV-C lamps are angled at 4 degrees to increase the dosage on the ceiling. This is because ceiling surfaces are harder to sanitize while the floor is easily mopped on a routine basis.
- Quartz glass is used for the UV lamps as it blocks less UV-C light than plain glass tubes, maximizing the actual output for the same wattage bulbs. Quartz glass also has increased strength to reduce chances of bulb breakage.
- Low pressure lamps are utilized since they produce virtually all of their output as UV-C light.
- UV output was designed to obtain greater than 99% reduction of typical viruses and bacteria in a 1-minute timeframe and on spores like C. diff in a 5-minute timeframe at a distance of 8 feet.
- Optional: A UV intensity sensor is available to both monitor the intensity and calculate dosage.

Economical:

The Torch is designed to be the lowest cost, high output UV generator available.

- Priced to allow users to purchase multiple Torches.



- Quartz glass is used for the UV lamps as it extends the bulb life by providing a better seal of the internal gasses.
- Low ozone production. The type of fused quartz used to make the body of the germicidal lamp determines the emission of the wavelength of the UV energy. Low ozone generating lamps transmit up to 90% of their energy at the 254nm wavelength and typically utilize a doped fused quartz that blocks the emission of 185nm energy.
- Solid state premium ballasts are used as they extend the bulb life by reducing the shock to the lamps when power is first turned on.
- Combined benefits of the UV lamps and ballasts extend the rated life of the lamps to over 16,000 hours.
- Replacement UV-C lamps are much less expensive than lamps from other manufacturers.
- A typical 15 minute exposure uses 0.16 kw-hours (kWh) of energy. At an average cost of 8 cents per kWh, the cost of a typical 15 minute exposure is 1.3 cents.
- No disposal of chemicals.
- No special lamp recycling required.

Easy to Operate:

- Easily operated with minimal training.
- No special room preparation is required.
- No chemicals to store and handle.
- Computer controlled for optimal performance.
- Motion sensors to detect motion in the room.
- Optional: UV Sensor to trend exposure data.
- Optional: Data is downloadable, enabling it to be exported for archival purposes.

Safe to Operate:

- Each Torch tower has an emergency stop button to inhibit a cycle or abort the process if pressed.
- The Torch UV Disinfection System must be manually reset if safety device is tripped. This prevents inadvertent restart of UV exposure as a further safety precaution.
- The Torch is started from a remote push button from outside of the room compared to other units which are started by pressing a button on the UV unit while still inside the room. This eliminates the risk of accidental exposure to UV for personnel operating the Torch, since the UV exposure is not initiated until all people are out of the room and the door is verified to be closed.



Torch Emergency Stop



- Four motion sensors are located on the tower to abort the UV exposure if motion is sensed in room.
-

Specifications:

- Torch Tower
68" H x 23" D x 23" W (1727mm H x 584mm D x 584mm W)
110-240 VAC, 6 Amps, 50/60 HZ
71 lbs (32 kg)
- Lamps are rated for 16,000 hours.
- Lamp type: 4-pin, low pressure, UVC Germicidal, low ozone
- Lamp quantity: 8
- Power cable: 15 feet, hospital grade
- Produces an intensity of approximately 12 mJ/cm² per minute (200 µw/cm²) at an 8-ft. distance.
- 3" diameter hospital grade wheels, resilient monprene.



Torch Tower

Design Features:

Protective Cover – A heavy duty cover is supplied with the Torch Tower to cover it when moving it around or storing it to better protect the lamps from damage.

Lamp Guard – A stainless steel protective lattice is incorporated into the Torch Tower to help protect the lamps from accidental breakage due to bumping hazards or items falling on it when the Protective Cover is not in place.

Options:

UV Sensor - An optional UV sensor is available to display and log dosage and archive data. Alarms and room numbers are logged as well as UV data if the UV Sensor option is chosen.

Bulbs:

ClorDiSys Solutions utilizes quartz lamps in the Torch UV System. Quartz is the premier material for UV producing lamps. ClorDiSys utilizes standard bulb lengths and ballasts. Our bulbs offer the best electrical efficiency by converting up to 40% of electrical power into to UV power. Our bulbs have a warm-up time of approx. 30 - 60 sec. With our LongLife+™ coating process, our low pressure mercury lamps have an operating life of up to 16,000 hours, maintaining an end-of-life UV-C output of 80%.



Used Bulb Waste Disposal

Our germicidal lamps are Toxicity Characteristic Leaching Procedure (TCLP) compliant. Lamps that PASS the TCLP test are considered as non-hazardous waste by the EPA.

In 1990 the EPA developed the TCLP test to simulate the effect of disposing waste in conventional landfills under complex environmental conditions. The method is designed to determine the mobility of toxic material in liquid, solid and multiphase waste. The EPA developed the TCLP to determine the toxicity of waste. The TCLP test does NOT measure the total mercury content but rather the potential of mercury to leach into groundwater if the waste is disposed of in a landfill. TCLP is designed to simulate the leaching that the waste will undergo if disposed of in a sanitary landfill. This test includes mercury, lead, cadmium, and other hazardous materials. Passing this test for mercury, for instance, requires a yield of less than 0.2 milligrams per liter upon completion of the test.

While lamps that pass TCLP may be classified as non-hazardous waste by the EPA, ClorDiSys Solutions and Clean Hospitals strongly encourage the recycling of spent germicidal lamps. Please contact your local environmental agency for assistance with lamp recycling or visit www.lamprecycle.org.

Appendix 1 – About UV-C

For the past 100 years science has recognized the bactericide effects of the ultraviolet area of the electromagnetic spectrum. Below are some key contributions over the years:

1855 Arloing and Daclaux demonstrated sunlight killed *Bacillus anthracis* and *Tyrophthrix scaber*

1877 Downes and Blunt reported bacteria were inactivated by sunlight – violet blue spectrum most effective

1889 Widmark confirmed UV rays from arc lamps were responsible for inactivation

1892 Geisler used a prism and heliostat to show sunlight and electric arc lamps are lethal to *Bacillus Typhosus*

1903 Banard and Morgan determined UV spectrum 226-328 nm is biocidal

1932 Ehris and Noethling isolated biocidal spectrum to 253.7 nm

1957 Riley proves effectiveness for Tb control

1994 CDC acknowledges UV effectiveness for Tb control

1999 WHO recommends UVGI for Tb control

2014 UV-C used as part of the terminal cleaning procedure within the Nebraska Biocontainment Unit upon ebola patient discharge

2020 UV-C Disinfection recommended for the disinfection of N95 masks and other PPE during SARS-CoV-2 pandemic.

The specific wavelengths responsible for the biocidal properties are situated between 240 - 280 nanometers (nm) with a peak wavelength at 265 nm. They are known as UV-C (see figure 1 & 2).




Fig. 1 - UV-C in the spectrum of electromagnetic radiation.




Fig. 2 - spectral energy distribution curve for germicidal action and spectral power distribution for low and medium pressure UV lamps.

UV-Action:

ClorDiSys' low-pressure, mercury-arc germicidal lamps are specially designed to produce the highest amounts of UV radiation - where 90% of energy is typically generated at 254nm. This radiation is very close to the peak of the germicidal effectiveness curve of 265nm, the most lethal wavelength to microorganisms. (See figure 2).

Our germicidal lamps are used extensively in the air purification markets and have been utilized in applications such as food and beverage, medical, HVAC (Heating, Ventilation and Air Conditioning), and pharmaceutical disinfection.

Our bulbs generate energy in the UV spectrum to destroy microorganisms: Microorganisms include several distinct groups of disease-causing germs, i.e. viruses, bacteria, fungi, algae and protozoa. The target of UV disinfection is the genetic material – nucleic acid. As UV light penetrates through the cell and is absorbed by the nucleic acids, a rearrangement of the genetic information occurs, interfering with the cells ability to reproduce. A cell that can't reproduce is considered dead; since it is unable to multiply to infectious numbers within a host. The maximum absorption of UV light by the nucleic acid, DNA, occurs at a wavelength of 265nm. The germicidal lamp emitting UV at 254nm is operating very close to the optimized wavelength the optimized wavelength for maximum absorption by nucleic acids.

Advantages of UV Radiation

Our process is environmentally friendly such that there are no dangerous or toxic chemicals that require specialized storage and/or handling and there are no concerns of overdosing. Since no chemicals are added to the air/water there are no process byproducts to be concerned with. Our equipment is cost effective with low initial capital cost and low operating costs. The process is effective since UV radiation offers immediate treatment process with no requirements for holding tanks or long retention/exposure times.



Safety

As UV-C provides radiation, it is not safe to be in the room while disinfection is taking place. UV-C is classified as “reasonably anticipated to be a human carcinogen” by the National Toxicology Program. It presents a hazard to skin and eyes, so direct exposure to UV-C is always to be avoided. UV-C is blocked by a number of materials, including glass (but not quartz glass) and most clear plastics, so it is possible to safely observe a UV-C system if you are looking through a window.

The process is environmentally friendly in that there are no dangerous or toxic chemicals that require specialized storage or handling. Since no chemicals are added to the air/water, there are no process byproducts to be concerned with. The UV bulbs do not require special handling or disposal either, making the system a green alternative to chemical disinfectants. UV-C provides residue free disinfection, so there is no concern over dangerous residues that need to be wiped down or neutralized after the disinfection occurs.

There has been concern with regard to the residual odors that have been noted after rooms are disinfected with ultraviolet light. Sometimes this smell is associated with ozone, a harmful gas. In reality, this odor is due to UV-C reacting with human dead skin cells and hair from dust in the room. Up to 80% of airborne dust in homes, offices, and other indoor environments is made up of dead human skin and hair. Skin and hair cells consist of keratin, a protein, while hair also contains cysteine, an amino acid. When high energy UV-C light hits keratin/cysteine molecules, it has enough power to break their internal chemical bonds creating smaller, sulfur-containing compounds that fall into the categories of thiols. The human nose is extremely sensitive to thiols and can detect them at concentrations as low as 1 part per billion. Concentrations of thiol molecules after a UV-C disinfection are negligible when compared to the published acceptable exposure limit. This means that any odor present after a UV-C disinfection is not dangerous, making the room immediately safe to enter after a UV-C disinfection has been performed.



Ultraviolet Dose

The degree of inactivation by ultraviolet radiation is directly related to the UV dose applied. The UV dose is the product of UV intensity [I] (expressed as energy per unit surface area) and exposure time [T]. Therefore: $DOSE = I \times T$

This dose is commonly expressed as millijoule per square centimeter (mJ/cm²).

The reduction of micro-organisms is classified using a logarithmic scale. A single log reduction is a 90% reduction of organisms. A two log reduction is a 99% reduction of organisms, followed by a three log reduction (99.9%), etc. The UV-C exposure dosage needed for each level of reduction is shown in the table along with the published reference where the data came from.

The Torch produces an intensity of approximately 12 mJ/cm² per minute (200 μw/cm²) at a 8-ft. distance.

UV Dose (mJ/cm²) for Various Reduction Levels

Spore	90%	99%	99.9%	99.99%	99.999%	99.9999%	Reference
Bacillus anthracis spores – Anthrax spores	24.32	48.64	72.96	97.28			UV-Light.co.UK
Bacillus magaterium sp. spores	2.73	5.46	8.19	10.92			UV-Light.co.UK
Bacillus subtilis ATCC6633(spores	36	48.6	61	78			Chang et al. 1985
Clostridioides difficile (C. diff) spores	6.0	12.0	18.0	24.0			UV-Light.co.UK
Bacterium							
Aeromonas salmonicida	1.5	2.7	3.1	5.9			Liltved and Landfald 1996
Aeromonas hydrophila ATCC7966	1.1	2.6	3.9	5	6.7	8.6	Wilson et al. 1992
Bacillus anthracis – Anthrax	4.52	9.04	13.56	18.08			UV-Light.co.UK
Bacillus magaterium sp. (veg.)	1.3	2.6	3.9	5.2			UV-Light.co.UK
Bacillus paratyphus	3.2	6.4	9.6	12.8			UV-Light.co.UK
Bacillus subtilis	5.8	11.6	17.4	23.2			UV-Light.co.UK
Campylobacter jejuni ATCC 43429	1.6	3.4	4	4.6	5.9		Wilson et al. 1992
Citrobacter diversus	5	7	9	11.5	13		Giese and Darby 2000
Citrobacter freundii	5	9	13				Giese and Darby 2000
Clostridium tetani	13.0	22.0					Light Sources Inc. 2014
Corynebacterium diphtheriae	3.37	6.74	10.11	13.48			UV-Light.co.UK
Ebertelia typhosa	2.14	4.28	6.42	8.56			UV-Light.co.UK
Escherichia coli O157:H7 CCUG 29193	3.5	4.7	5.5	7			Sommer et al. 2000
Escherichia coli O157:H7	<2	<2	2.5	4	8	17	Yaun et al. 2003
Halobacterium elongate ATCC33173	0.4	0.7	1				Martin et al. 2000
Halobacterium salinarum ATCC43214	12	15	17.5	20			Martin et al. 2000
Klebsiella pneumoniae	12	15	17.5	20			Giese and Darby 2000
Klebsiella terrigena ATCC33257	4.6	6.7	8.9	11			Wilson et al. 1992

UV Dose (mJ/cm²) for Various Reduction Levels

Spore	90%	99%	99.9%	99.99%	99.999%	99.9999%	Reference
<i>Legionella pneumophila</i> ATCC33152	1.9	3.8	5.8	7.7	9.6		Oguma et al. 2004
<i>Leptospira</i> canicola – infectious Jaundice	3.15	6.3	9.45	12.6			UV-Light.co.UK
<i>Micrococcus candidus</i>	6.05	12.1	18.15	24.2			UV-Light.co.UK
<i>Micrococcus sphaeroides</i>	1.0	2.0	3.0	4.0			UV-Light.co.UK
<i>Mycobacterium tuberculosis</i>	6.2	12.4	18.6	24.8			UV-Light.co.UK
MRSA	3.2	6.4	9.6	12.8			UV-Light.co.UK
<i>Neisseria catarrhalis</i>	4.4	8.8	13.2	17.6			UV-Light.co.UK
<i>Phytomonas tumefaciens</i>	4.4	8.8	13.2	17.6			UV-Light.co.UK
<i>Proteus vulgaris</i>	3.0	6.0	9.0	12.0			UV-Light.co.UK
<i>Pseudomonas stutzeri</i>	100	150	195	230			Joux et al. 1999
<i>Pseudomonas aeruginosa</i>	5.5	11.0	16.5	22.0			UV-Light.co.UK
<i>Pseudomonas fluorescens</i>	3.5	7.0	10.5	14.0			UV-Light.co.UK
<i>Salmonella anatum</i> (from human feces)	7.5	12	15				Tosa and Hirata 1998
<i>Salmonella derby</i> (from human feces)	3.5	7.5					Tosa and Hirata 1998
<i>Salmonella enteritidis</i>	4.0	8.0	12.0	16.0			UV-Light.co.UK
<i>Salmonella infantis</i> (from human feces)	2	4	6				Tosa and Hirata 1998
<i>Salmonella paratyphi</i> – Enteric fever	3.2	6.4	9.6	12.8			UV-Light.co.UK
<i>Salmonella typhosa</i> – Typhoid fever	2.15	4.3	6.45	8.6			UV-Light.co.UK
<i>Salmonella typhimurium</i>	8.0	16.0	24.0	32.0			UV-Light.co.UK
<i>Sarcina lutea</i>	19.7	39.4	59.1	78.8			UV-Light.co.UK
<i>Serratia marcescens</i>	2.42	4.84	7.26	9.68			UV-Light.co.UK
<i>Shigella dysenteriae</i> – Dysentery	2.2	4.4	6.6	8.8			UV-Light.co.UK
<i>Shigella flexneri</i> – Dysentery	1.7	3.4	5.1	6.8			UV-Light.co.UK
<i>Shigella paradysenteriae</i>	1.68	3.3	5.04	6.72			UV-Light.co.UK
<i>Shigella sonnei</i> ATCC9290	3.2	4.9	6.5	8.2			Chang et al. 1985
<i>Spirillum rubrum</i>	4.4	8.8	13.2	17.6			UV-Light.co.UK
<i>Staphylococcus albus</i>	1.84	3.68	5.52	7.36			UV-Light.co.UK
<i>Staphylococcus aureus</i>	2.6	5.2	7.8	10.4			UV-Light.co.UK
<i>Staphylococcus hemolyticus</i>	2.16	4.32	6.48	8.64			UV-Light.co.UK
<i>Staphylococcus lactis</i>	6.15	12.3	18.45	24.6			UV-Light.co.UK
<i>Streptococcus faecalis</i> ATCC29212	6.6	8.8	9.9	11.2			Chang et al. 1985
<i>Streptococcus viridans</i>	2.0	4.0	6.0	8.0			UV-Light.co.UK
<i>Vibrio anguillarum</i>	0.5	1.2	1.5	2			Liltved and Landfald 1996
<i>Vibrio comma</i> – Cholera	3.375	6.75	10.125	13.5			UV-Light.co.UK
<i>Vibrio natriegens</i>	37.5	75	100	130	150		Joux et al. 1999
<i>Yersinia enterocolitica</i> ATCC27729	1.7	2.8	3.7	4.6			Wilson et al. 1992
<i>Yersinia ruckeri</i>	1	2	3	5			Liltved and Landfald 1996
Yeasts							
Brewers yeast	3.3	6.6	9.9	13.2			UV-Light.co.UK
Common yeast cake	6.0	12.0	18.0	24.0			UV-Light.co.UK
<i>Saccharomyces cerevisiae</i>	6.0	12.0	18.0	24.0			UV-Light.co.UK



UV Dose (mJ/cm²) for Various Reduction Levels

Spore	90%	99%	99.9%	99.99%	99.999%	99.9999%	Reference
<i>Saccharomyces ellipsoideus</i>	6.0	12.0	18.0	24.0			UV-Light.co.UK
<i>Saccharomyces</i> spores	8.0	16.0	24.0	32.0			UV-Light.co.UK
Molds							
<i>Aspergillus flavus</i>	60.0	120.0	180.0	240.0			UV-Light.co.UK
<i>Aspergillus glaucus</i>	44.0	88.0	132.0	176.0			UV-Light.co.UK
<i>Aspergillus niger</i>	132.0	264.0	396.0	528.0			UV-Light.co.UK
<i>Mucor racemosus</i> A	17.0	34.0	51.0	68.0			UV-Light.co.UK
<i>Mucor racemosus</i> B	17.0	34.0	51.0	68.0			UV-Light.co.UK
<i>Oospora lactis</i>	5.0	10.0	15.0	20.0			UV-Light.co.UK
<i>Penicillium digitatum</i>	44.0	88.0	132.0	176.0			UV-Light.co.UK
<i>Penicillium expansum</i>	13.0	26.0	39.0	52.0			UV-Light.co.UK
<i>Penicillium roqueforti</i>	13.0	26.0	39.0	52.0			UV-Light.co.UK
<i>Rhizopus nigricans</i>	111.0	222.0	333.0	444.0			UV-Light.co.UK
Protozoan							
<i>Chlorella Vulgaris</i>	13.0	26.0	39.0	52.0			UV-Light.co.UK
<i>Cryptosporidium hominis</i>	3	5.8					Johnson et al. 2005
<i>Cryptosporidium parvum</i>	2.4	<5	5.2	9.5			Craik et al. 2001
<i>Cryptosporidium parvum</i> , oocysts, tissue culture assay	1.3	2.3	3.2				Shin et al. 2000
<i>Encephalitozoon cuniculi</i> , microsporidia	4	9	13				Marshall et al. 2003
<i>Encephalitozoon hellem</i> , microsporidia	8	12	18				Marshall et al. 2003
<i>Encephalitozoon intestinalis</i> , microsporidia	<3	3	<6	6			Huffman et al. 2002
<i>Giardia lamblia</i>	<10	~10	<20				Campbell et al. 2002
<i>Giardia muris</i>	<10	<10	<25	~60			Belosevic et al. 2001
Nematode Eggs	45.0	90.0	135.0	180.0			UV-Light.co.UK
Paramecium	11.0	22.0	33.0	44.0			UV-Light.co.UK



The following table shows the reduction values for various viruses.

UV Dose (mJ/cm ²) for Various Reduction Levels								
Virus	Host	90%	99%	99.9%	99.99%	99.999%	99.9999%	Reference
Adenovirus type 15	A549 cell line (ATCC CCL-	40	80	122	165	210		Thompson et al. 2003
Adenovirus type 2	PLC / PRF / 5	40	78	119	160	195	235	Gerba et al. 2002
B40-8 (Phage)	B. Fragilis	11	17	23	29	35	41	Sommer et al. 2001
Bacteriophage – E. Coli		2.6	5.2	7.8	104.0			UV-Light.co.UK
Calicivirus canine	MDCK cell line	7	15	22	30	36		Husman et al. 2004
Calicivirus feline	CRFK cell line	5	15	23	30	39		Thurston-Enriquez et al. 2003
Coxsackievirus B3	BGM cell line	8	16	24.5	32.5			Gerba et al. 2002
Coxsackievirus B5	BGM cell line	9.5	18	27	36			Gerba et al. 2002
Echovirus I	BGM cell line	8	16.5	25	33			Gerba et al. 2002
Echovirus II	BGM cell line	7	14	20.5	28			Gerba et al. 2002
Hepatitis A HM175	FRhK-4 cell	5.1	13.7	22	29.6			Wilson et al. 1992
Infectious Hepatitis		5.8	11.6	17.4	232.0			UV-Light.co.UK
Influenza		3.4	6.8	10.2	136.0			UV-Light.co.UK
MS2 (Phage)	E. coli		45	75	100	125	155	Thompson et al. 2003
Norovirus		10	16	22	26	30		Lee et al. 2008
Parvovirus		2.2	4.6					Cornelis et al. 1982
PHI X 174 (Phage)	E. coli WG 5	3	5	7.5	10	12.5	15	Sommer et al. 2001
Poliovirus – Poliomyelitis		3.15	6.3	9.45	126.0			UV-Light.co.UK
Poliovirus 1	CaCo2 cell-line (ATCC HTB37)	7	17	28	37			Thompson et al. 2003
PRD-1 (Phage)	S. typhimurium	9.9	17.2	23.5	30.1			Meng and Gerba 1996
Reovirus Type 1 Lang strain	N/A	16	36					Harris et al. 1987
Reovirus-3	Mouse L-60	11.2	22.4					Rauth 1965
Rotavirus	MA104 cells	20	80	140	200			Caballero et al. 2004
Rotavirus SA-11	MA-104 cell	9.1	19	26	36	48		Wilson et al. 1992
SARS-CoV-2	N/A		5				22	Boston University, 2020
Staphylococcus aureus phage A	Staphylococcus aureus 994	8	17	25	36	47		Sommer et al. 1989
Tobacco mosaic	N/A	240.0	440.0					Light Sources Inc. 2014



Appendix 2 – Persistence of Bacteria

(As compiled via a Google Search)

Persistence of Clinically Relevant Bacteria on Dry Inanimate Surfaces ¹	
Organism	Persistence
Acinetobacter spp.	3 days - 5 months
Bordetella pertussis	3-5 days
Campylobacter jejuni	Up to 6 days
Clostridium difficile (spores)	5 months
Chlamydia pneumoniae	Up to 30 hours
Chlamydia psittaci	15 days
Corynebacterium diphtheria	7 days – 6 months
Corynebacterium pseudotuberculosis	1-8 days
Escherichia coli	1.5 hours – 16 months
Enterococcus spp. including VRE and VSE	5 days – 4 months
Haemophilus influenza	12 days
Helicobacter pylori	Up to 90 minutes
Klebsiella spp.	2 hours – 30 months
Listeria spp.	1 day – 4 months
Mycobacterium bovis	Up to 2 months
Mycobacterium tuberculosis	1 day – 4 months
Neisseria gonorrhoeae	1-3 days
Proteus vulgaris	1-2 days
Pseudomonas aeruginosa	6 hours – 16 months; 5 weeks on dry floor
Salmonella typhi	6 hours – 4 weeks
Salmonella typhimurium	10 days – 4.2 years
Salmonella spp.	1 day
Serratia marcescens	3 days – 2 months; 5 weeks on dry floor
Shigella spp.	2 days – 5 months
Staphylococcus aureus, including MRSA	7 days – 7 months
Streptococcus pneumoniae	1-20 days
Streptococcus pyogenes	3 days – 6.5 months
Vibrio cholera	1-7 days

References:

- Amoah, K., Craik, S., Smith, D.W. and Belosevic, M. 2005. Inactivation of *Cryptosporidium* oocysts and *Giardia* cysts by ultraviolet light in the presence of natural particulate matter, *AQUA, J. Wat. Supply* 54(3): 165-178.
- Ballester, N.A. and Malley, J.P. 2004. Sequential disinfection of adenovirus type 2 with UV-chlorinechloramine, *J. Amer. Wat. Works Assoc.*, 96(10): 97-102.
- Batch, L.F., Schulz, C.R. and Linden, K.G. 2004. Evaluating water quality effects on UV disinfection of MS2 coliphage, *J. Amer. Wat. Works Assoc.*, 96(7): 75-87.
- Battigelli, D.A., Sobsey, M.D. and Lobe, D.C. 1993. The inactivation of Hepatitis A virus and other model viruses by UV irradiation, *Wat. Sci. Tech.*, 27(3-4): 339-342.
- Belosevic, M., Craik, S.A., Stafford, J.L. Neumann, N.E., Kruithof, J. and Smith, D.W. 2001. Studies on the resistance/reaction of *Giardia muris* cysts and *C. parvum* oocysts exposed to medium-pressure ultraviolet radiation, *FEMS Microbiol. Lett.*, 204(1): 197-204.
- Bolton J.R. and Linden, K.G. 2003. Standardization of methods for fluence (UV Dose) determination in benchscale UV experiments. *J. Environ. Eng.* 129(3): 209-216.
- Bukhari, Z., Abrams, F. and LeChevallier, M. 2004. Using ultraviolet light for disinfection of finished water, *Water Sci. Tech.*, 50(1): 173-178.
- Caballero, S., Abad, F.X., Loisy, F., Le Guyader, F.S., Cohen, J., Pinto, R.M. and Bosch, A. 2004. Rotavirus virus-like particles as surrogates in environmental persistence and inactivation studies, *Appl. Env. Microbiol.* 70(7): 3904-3909.
- Campbell, A.T. and Wallis, P. 2002. The effect of UV irradiation on human-derived *Giardia lamblia* cysts, *Wat. Res.*, 36(4): 963- 969.
- Carlson, D.A., Seabloom, R.W., DeWalle, F.B., Wetzler, T.F., Engeset, J., Butler, R., Wangsuphachart, S. and Wang, S. 1985. Ultraviolet disinfection of water for small water supplies. US EPA Report No. EPA/600/S2-85/092.
- Cass AL, Kelly JW, Probst JC, Addy CL, McKeown RE. Identification of device-associated infections utilizing administrative data. *American Journal of Infection Control*. 2013; Published online 17 June 2013.
- Chang, J.C.H., Osoff, S.F., Lobe, D.C., Dorfman, M.H., Dumais, C.M., Qualls, R.G. and Johnson, J.D. 1985. UV inactivation of pathogenic and indicator microorganisms, *Appl. Environ. Microbiol.*, 49(6): 1361-1365.
- Clancy, J.L., Bukhari, Z., Hargy, T.M., Bolton, J.R., Dussert, B.W. and Marshall, M.M. 2000. Using UV to inactivate *Cryptosporidium* – Even extremely low dosages of ultraviolet light can be highly effective for inactivating *Cryptosporidium* oocysts, *J. Amer. Wat. Works Assoc.*, 92(9): 97-104.
- Clancy, J.L., Marshall, M.M., Hargy, T.M. and Korich, D.G. 2004. Susceptibility of five strains of *Cryptosporidium parvum* oocysts to UV light, *J. Amer. Wat. Works Assoc.*, 96(3), 84-93.
- Cornelis, J.J., Su, Z.Z., and Rommelaere, J. 1982. Direct and Indirect Effects of Ultraviolet Light on the Mutagenesis of Parvovirus H-1 in Human Cells, *The EMBO Journal*, 1(6):693-699.
- Craik, S.A., Finch, G.R., Bolton, J.R. and Belosevic, M. 2000. Inactivation of *Giardia muris* cysts using medium- pressure ultraviolet radiation in filtered water, *Wat. Res.*,34(18): 4325-4332.
- Craik, S.A., Weldon, D., Finch, G.R., Bolton, J.R. and Belosevic, M. 2001. Inactivation of *Cryptosporidium parvum* oocysts using medium- and low-pressure ultraviolet radiation, *Wat. Res.*, 35(6): 1387-1398.
- Faires MC, Pearl DL, Berke O, Reid-Smith RJ, Weese JS. The identification and epidemiology of methicillin-resistant *Staphylococcus aureus* and *Clostridium difficile* in patient rooms and the ward environment. *BMC Infectious Diseases* 2013, 13:342
- Gerba, C.P., Gramos, D.M. and Nwachuku, N. 2002. Comparative inactivation of enteroviruses and adenovirus 2 by UV light, *Appl. Environ. Microbiol.*,68(10): 5167-5169.
- Giese, N. and Darby, J. 2000. Sensitivity of microorganisms to different wavelengths of UV light: implications on modeling of medium pressure UV systems, *Wat. Res.*, 34(16): 4007-4013.
- Grohskopf LA, Sinkowitz-Cochran RL, Garrett Dom et al. A national point-prevalence survey of pediatric intensive care unit-acquired infections in the United States. *Journal of Pediatrics*. 2002; 140, 432-438.
- Harris, G.D., Adams, V.D., Sorensen, D.L. and Curtis, M.S. 1987. Ultraviolet inactivation of selected bacteria and viruses with photoreactivation of the bacteria, *Wat. Res.*,21(6): 687-692.
- Hayes, S.L., Rice, E.W., Ware, M.W. and Schaefer III, F.W.2003. Low pressure ultraviolet studies for inactivation of *Giardia muris* cysts, *J. Appl. Microbiol.*, 94(1): 54-59.
- Hijnen, W.A.M., Beerendonk, E.F. and Medema, G.J. 2006. Inactivation credit of UV radiation for viruses, bacteria and protozoan (oo)cysts in water; a review, *Wat. Res.*,40(1): 3-22.
- Hoyer, O. 1998. Testing performance and monitoring of UV systems for drinking water disinfection, *Wat. Supply*,16(1-2): 424-429.
- Huffman, D.E., Gennaccaro, A., Rose, J.B. and Dussert, B.W. 2002. Low- and medium-pressure UV inactivation of microsporidia *Encephalitozoon intestinalis*, *Wat. Res.*,36(12): 3161-3164.
- Husman, A.M.D., Bijkerk, P., Lodder, W., Van den Berg, H., Pribil, W., Cabaj, A., Gehringer, P., Sommer, R. and Duizer, E. 2004. Calicivirus inactivation by nonionizing 253.7-nanometer-wavelength (UV) and ionizing (Gamma) radiation, *Appl. Environ. Microbiol.*, 70(9):5089-5093.
- Infection Control Breaches in the Operating Room. *Infection Control Today*. 2013.
- Johnson, A.M., Linden, K., Ciociola, K.M., De Leon, R., Widmer, G. and Rochelle, P.A. 2005. UV inactivation of *Cryptosporidium hominis* as measured in cell culture, *Appl. Environ. Microbiol.*, 71(5): 2800-2802.
- Joux, F., Jeffrey, W.H., Lebaron, P. and Mitchell, D. L. 1999. Marine bacterial isolates display diverse responses to UV-B radiation, *Appl. Environ. Microbiol.*, 65(9):3820-3827.

- Karanis, P., Maier, W.A., Seitz, H.M. and Schoenen, D. 1992. UV sensitivity of protozoan parasites, *Aqua*, 41:95-100.
- Lazarova, V. and Savoye, P. 2004. Technical and sanitary aspect of wastewater disinfection by ultraviolet irradiation for landscape irrigation, *Wat. Sci. Technol.*, 50(2): 203-209.
- Lee, J., Zoh, K., and Ko, G. 2008. Inactivation and UV Disinfection of Murine Norovirus with TiO₂ under Various Environmental Conditions, *Applied and Environmental Microbiology.*, 74(7):2111-2117.
- Light Sources Inc and American Ultraviolet Company. UV Irradiation Dosage Table. Accessed from <http://www.americanairandwater.com/uv-facts/uv-dosage.htm>. Accessed on 3-26-2014
- Liltved, H. and Landfald, B. 1996. Influence of liquid holding recovery and photoreactivation on survival of ultraviolet-irradiated fish pathogenic bacteria, *Wat. Res.*, 30(5): 1109-1114.
- Linden, K.G., Batch, L. and Schulz, C. 2002a. UV disinfection of filtered water supplies: water quality impacts on MS2 dose-response curves, *Proceedings Amer. Wat. Works Assoc. Annu. Conf.*, Amer. Wat. Works Assoc., Denver, CO.
- Linden, K.G., Shin, G.-A., Faubert, G., Cairns, W. and Sobsey, M.D. 2002b. UV disinfection of *Giardia lamblia* cysts in water, *Environ. Sci. Technol.*, 36(11): 2519-2522.
- Lucado, J., Paez, K., Andrews, R., Steiner, C.. Adult Hospital Stays with Infections Due to Medical Care, 2007.HCUP Statistical Brief #94. August 2010. Agency for Healthcare Research and Quality, Rockville, MD.
- Mamane-Gravetz, H. and Linden, K.G. 2004. UV disinfection of indigenous aerobic spores: Implications for UV reactor validation in unfiltered waters, *Wat. Res.*, 38(12): 2898-2906.
- Marshall, M.M., Hayes, S., Moffett, J., Sterling, C.R. and Nicholson, W.L. 2003. Comparison of UV inactivation of three *Encephalitozoon* species with that of spores of two DNA repair-deficient *Bacillus subtilis* biosimetry strains, *Appl. Environ. Microbiol.*, 69(1): 683-685.
- Martin, E.L., Reinhardt, R.L., Baum, L.L., Becker, M.R., Shaffer, J.J. and Kokjohn, T.A. 2000. The effects of ultraviolet radiation on the moderate halophile *Halomonas elongata* and the extreme halophile *Halobacterium salinarum*, *Can. J. Microbiol.*, 46(2): 180-187.
- Maya, C., Beltran, N., Jimenez, B. and Bonilla, P. 2003. Evaluation of the UV disinfection process in bacteria and amphizoic amoebae inactivation, *Wat. Sci. Technol.: Wat. Supply*, 3(4): 285-291.
- Meng, Q.S. and Gerba, C.P. 1996. Comparative inactivation of enteric adenoviruses, poliovirus and coliphages by ultraviolet irradiation, *Wat. Res.*, 30(11):2665-2668.
- Mofidi, A.A., Meyer, E.A., Wallis, P.M., Chou, C.I., Meyer, B.P., Ramalingam, S. and Coffey, B.M. 2002. The effect of UV light on the inactivation of *Giardia lamblia* and *Giardia muris* cysts as determined by animal infectivity assay, *Wat. Res.*, 36(8): 2098-2108.
- Morita, S., Namikoshi, A., Hirata, T., Oguma, K., Katayama, H., Ohgaki, S., Motoyama, N. and Fujiwara, M. 2002. Efficacy of UV irradiation in inactivating *C. parvum* oocysts, *Appl. Environ. Microbiol.*, 68(11):5387-5393.
- Nieuwstad, T.J. and Havelaar, A.H. 1994. The kinetics of batch ultraviolet inactivation of bacteriophage MS2 and microbiological calibration of an ultraviolet pilot plant, *J. Environ. Sci. Health*, A29(9): 1993-2007.
- Oguma, K., Katayama, H. and Ohgaki, S. 2002. Photoreactivation of *E. coli* after low- or medium- pressure UV disinfection determined by an endonuclease sensitive site assay, *Appl. Environ. Microbiol.*, 68(12), 6029-6035.
- Oguma, K., Katayama, H. and Ohgaki, S. 2004. Photoreactivation of *Legionella pneumophila* after inactivation by low- or medium-pressure ultraviolet lamp, *Wat. Res.*, 38(11): 2757-2763.
- Oppenheimer, J.A., Hoagland, J.E., Laine, J.-M., Jacangelo, J.G. and Bhamrah, A. 1993. Microbial inactivation and characterization of toxicity and by-products occurring in reclaimed wastewater disinfected with UV radiation, *Conf. on Planning, Design and Operation of Effluent Disinfection Systems*, Whippany, NJ, May 23-25, 1993,
- Pyrek, K. Pathogen Persistence, Transmission and Cross-Contamination Prevention. VIRGO Publishing. Aug 2012.
- Wat. Environ. Fed., Alexandria, VA Otaki, M., Okuda, A., Tajima, K., Iwasaki, T., Kinoshita, S. and Ohgaki, S. 2003. Inactivation differences of microorganisms by low pressure UV and pulsed xenon lamps, *Wat. Sci. Technol.*, 47(3): 185-190.
- Rauth, A.M. 1965. The physical state of viral nucleic acid and the sensitivity of viruses to ultraviolet light, *Biophys. J.*, 5: 257-273.
- Rice, E.W. and Hoff, J.C. 1981. Inactivation of *Giardia lamblia* cysts by ultraviolet irradiation, *Appl. Environ. Microbiol.*, 42(3): 546-547.
- Shin, G.-A., Linden, K.G. and Sobsey, M.D. 2000. Comparative inactivation of *Cryptosporidium parvum* oocysts and coliphage MS2 by monochromatic UV radiation, *Proceedings of Disinfection 2000: Disinfection of Wastes in the New Millennium*, New Orleans, Water Environment Federation, Alexandria, VA.
- Shin, G.-A., Linden, K.G., Arrowood, M.J. and Sobsey, M.D. 2001. Low-pressure UV inactivation and DNA repair potential of *C. parvum* oocysts, *Appl. Environ. Microbiol.*, 67(7): 3029-3032.
- Shin, G.A., Linden, K.G. and Sobsey, M.D. 2005. Low pressure ultraviolet inactivation of pathogenic enteric viruses and bacteriophages, *J. Environ. Engr. Sci.*, 4: S7-S11.
- Sommer, R., Weber, G., Cabaj, A., Wekerle, J., Keck, G., and Schauburger, G. 1989. UV inactivation of micro-organisms in water. *Zbl. Hyg.* 189: 214-224.
- Sommer, R., Haider, T., Cabaj, A., Pribil, W. and Lhotsky, M. 1998. Time dose reciprocity in UV disinfection of water, *Water Sci. Technol.*, 38(12): 145-150.
- Sommer, R., Cabaj, A., Sandu, T. and Lhotsky, M. 1999. Measurement of UV radiation using suspensions of microorganisms, *J. Photochem. Photobiol.*, 53(1-3): 1-5.
- Sommer, R., Lhotsky, M., Haider, T. and Cabaj, A. 2000. UV inactivation, liquid-holding recovery, and photoreactivation of *E. coli* O157 and other pathogenic *E. coli* strains in water, *J. Food Protection*, 63(8): 1015-1020.
- Sommer, R., Pribil, W., Appelt, S., Gehringer, P., Eschweiler, H., Leth, H., Cabaj, A. and Haider, T. 2001. Inactivation of bacteriophages in water by means of non-ionizing (UV-253.7



- nm) and ionizing (gamma) radiation: A comparative approach, *Wat. Res.*, 35(13): 3109- 3116.
- Srinivasan MD, Arjun, American Recovery and Reinvestment Act Epidemiology and Laboratory Capacity (ELC) for Infectious Disease Program Healthcare-Associated Infections (HAIs) Grantee Meeting CDR Oct 19-20 2009.
- Thurston-Enriquez, J.A. , Haas, C.N. , Jacangelo, J. , Riley, K. and Gerba, C.P. 2003. Inactivation of feline calicivirus and adenovirus type 40 by UV radiation, *Appl. Environ. Microbiol.*, 69(1): 577-582.
- Thompson, S.S., Jackson, J.L., Suva-Castillo, M., Yanko, W.A., Jack, Z.E., Kuo, J., Chen, C.L., Williams, F.P. and Schnurr, D.P. 2003. Detection of infectious human adenoviruses in tertiary-treated and ultraviolet-disinfected wastewater, *Wat. Environ. Res.*, 75(2): 163-170.
- Tosa, K. and Hirata, T. 1998. HRWM-39: Photoreactivation of Salmonella following UV disinfection, IAWQ 19th Biennial International Conference, Vol. 10, Health- Related Water Microbiology.
- Tosa, K. and Hirata, T. 1999. Photoreactivation of enterohemorrhagic *E. coli* following UV disinfection, *Wat. Res.*, 33(2): 361-366.
- Tree, J.A., Adams, M.R. and Lees, D.N. 1997. Virus inactivation during disinfection of wastewater by chlorination and UV irradiation and the efficacy of F+ bacteriophage as a 'viral indicator', *Wat. Sci. Technol.*, 35(11-12): 227-232.
- Tree, J.A., Adams, M.R. and Lees, D.N. 2005. Disinfection of feline calicivirus (a surrogate for Norovirus) in wastewaters, *J. Appl. Microbiol.*, 98: 155-162.
- UV-Light.co.UK, UV Light Technology Limited, <https://www.uv-light.co.uk/uv-dose-required-for-inactivation-of-viruses-bacteria-moulds-etc/> accessed on 2-20-2018.
- Wiedenmann, A. , Fischer, B., Straub, U., Wang, C.-H., Flehmig, B. and Schoenen, D. 1993. Disinfection of Hepatitis A virus and MS-2 coliphage in water by ultraviolet irradiation: Comparison of UV-susceptibility, *Wat. Sci. Tech.*, 27(3-4): 335-338.
- Wiener-Well Y, Galuty M, Rudensky B, Schlesinger Y, Attias D, Yinnon AM. Nursing and physician attire as possible source of nosocomial infections. *American Journal of Infection Control*. 2011 Sep;39(7):555-9.
- Wilson, B.R., Roessler, P.F., Van Dellen, E., Abbaszadegan, M. and Gerba, C.P. 1992. Coliphage MS-2 as a UV water disinfection efficacy test surrogate for bacterial and viral pathogens, Proceedings, Water Quality Technology Conference, Nov 15-19, 1992, Toronto, Canada, pp. 219-235, Amer. Wat. Works Assoc., Denver, CO.
- Wu, Y., Clevenger, T. and Deng, B. 2005. Impacts of goethite particles on UV disinfection of drinking water, *Appl. Environ. Microbiol.*, 71(7): 4140-4143.
- Yaun, B.R., Sumner, S.S., Eifert, J.D. and Marcy, J.E. 2003. Response of *Salmonella* and *E. coli* O157:H7 to UV energy, *J. Food Protection*, 66(6): 1071-1073.
- Zimmer, J.L. and Slawson, R.M. 2002. Potential repair of *E. coli* DNA following exposure to UV radiation from both medium- and low-pressure UV sources used in drinking water treatment, *Appl. Environ. Microbiol.*, 68(7): 3293-3299.
- Zimmer, J.L., Slawson, R.M. and Huck, P.M. 2003. Inactivation and potential repair of *C. parvum* following low- and medium-pressure ultraviolet irradiation, *Wat. Res.*, 37(14): 3517-3523.