



"The Chlorine Dioxide People"

Providing you with gaseous chlorine dioxide solutions for your decontamination needs

Portable Decontamination Cart

Description:

The ClorDiSys Portable Decontamination Cart offers the ability to decontaminate small equipment and components inside a small chamber rather than in an entire room or large chamber. This allows for a quicker, more cost effective cycle for decontaminating a few items by limiting the amount of consumables used. Components sterilized within the chamber can then be transported throughout the facility in a sterile environment.

With both the Chlorine Dioxide Generator and the Decontamination Cart being easily portable, a complete decontamination station can be setup anywhere in your facility. With this system, you can decontaminate infected or contaminated items without taking them out of the room, minimizing risk of spreading the contamination. It provides a rapid and highly effective method to sterilize computers, electronics, medical devices, sterile products, instruments, animal cages, tools, supplies and other components at ambient temperatures.

Constructed out of polypropylene, the PDC comes equipped with all necessary connections to interface with our Chlorine Dioxide Gas generators, including gas inject and gas concentration sampling ports, and aeration ports as well. It also allows for the gas decontamination of the HEPA filters from caging racks which cannot go into an autoclave and much more. The Portable Decontamination Cart can be made in a variety of sizes, depending on your facility's needs.



Custom sizing is available

Decontamination time can be under **1 hour** for a 20 ft³ (0.566 m³) chamber

Costs of under **\$5.00 in consumables** per cycle

Call us for more information about the power of a TRUE GAS
908.236.4100



P.O. Box 549, Lebanon, NJ 08833-0549 Tel:(908) 236-4100 Fax:(908) 236-2222
www.clordisys.com

44-040 Rev 10-2014