



"The Chlorine Dioxide People"

Providing you with gaseous chlorine dioxide solutions for your decontamination needs

Portable Decontamination Chamber

The Portable Chlorine Dioxide Gas Decontamination Chamber is an easily transportable, sealed chamber designed for use in any governmental, pharmaceutical, manufacturing, laboratory, research or surgical setting. It is used in conjunction with a CD Generator to provide a rapid and highly effective method to sterilize computers, electronics, medical devices, sterile products, instruments, and components and transport them under aseptic conditions to a different part of your facility. It also offers a way to sterilize components within the room where they are kept to minimize the chance for cross-contamination. The Portable Decon Chamber has HEPA filtered air intakes to keep components sterile after decontamination.

The Portable Decon Chamber contains removable shelving inside to handle multiple or large components. It can connect to your facility's exhaust system for direct aeration or utilize an integrated carbon scrubber to remove the gas from the chamber.



*Dimensions: 60" L x 60" H x 30" D
(152 cm x 152 cm x 76 cm)*

Custom sizing is available to fit your facility's needs

Decontamination time can be under **1.5 hours** for a 63 ft³ (1.8 m³) chamber

Costs of under **\$7 in consumables** per cycle

Call us for more information about the power of a TRUE GAS
908.236.4100



P.O. Box 549, Lebanon, NJ 08833-0549 Tel:(908) 236-4100 Fax:(908) 236-2222
www.clordisys.com

14-048 Rev 10-2014